

Python On Resonance (PyOR)

Everybody can simulate NMR

Author: Vineeth Thalakkotloor

Email: vineethfrancis.physics@gmail.com

Tutorial 5: Evolution of Density Matrix in Hilbert Space Part 1

In previous tutorial, we saw how to generate initial density matrix for two spin system a electron and a proton at thermal equilibrium. In this tutorial we try to evolve the density matrix in time by solving Liouville-Von Neumann equation in Hilbert Space for a single spin half system

Load Python packages and define path to the source file "PythonOnResonance.py"

```
In [67]: pathSource = '/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source'
```

```
In [68]: from IPython.display import display, HTML
display(HTML("<style>.container { width:100% !important; }</style>"))
import sys
sys.path.append(pathSource)

import PythonOnResonance as PyOR

import time
import numpy as np
import matplotlib.pyplot as plt
from matplotlib import rc
%matplotlib notebook
import sympy as sp
from sympy import *
```

Generating Spin System

```
In [69]: """
Define Spin quantum numbers of your spins in "Slist1".
Slist1[0] is spin of first particle and Slist1[1] is spin of second particle.
""";

Slist1 = [1/2]
```

```
In [70]: """
Define Planck constant equals 1.
Because NMR spectroscopists are more interested to write Energy in frequency units.
if False then hbarEQ1 = hbar
""";

hbarEQ1 = True
```

```
In [71]: """
Generate Spin Operators
""";

System = PyOR.Numerical_MR(Slist1,hbarEQ1)

"""
Sx, Sy and Sz Operators
""";
Sx,Sy,Sz = System.SpinOperator()

"""
S+ and S- Operators
""";
Sp,Sm = System.PMoperators(Sx,Sy)
```

Zeeman Hamiltonian in Lab Frame

```
In [72]: """
Gyromagnetic Ratio
Gamma = [Gyromagnetic Ratio spin 1, Gyromagnetic Ratio spin 1, ...]
""";
Gamma = [System.gammaH1]

"""
Define the field of the spectrometer, B0 in Tesla.
""";
B0 = 9.4

"""
Define the chemical Shift of individual spins
Offset = [chemical Shift spin 1, chemical Shift spin 1, ..]
""";
Offset = [20] # Offset frequency in Hz

"""
Function "LarmorF" give the list Larmor frequencies of individual spins in lab frame
""";
LarmorF = System.LarmorFrequency(Gamma,B0,Offset)

Hz = System.Zeeman(LarmorF,Sz)
```

Larmor Frequency in MHz: [-400.22803765]

Initialize Density Matrix

```
In [73]: """
We will generate Initial Density Matrix in two ways:
First we will generate a density matrix as we prefer say, Sz.
Second we will create density matrix at thermal equilibrium

First Case
""";

Thermal_DensMatrix = False

if Thermal_DensMatrix:
    Hz_EnUnit = System.Convert_FreqUnitsTOEnergy(Hz)
    HT_approx = False # High Temperature Approximation is False
    T = 300 # Temperature in Kelvin
```

```

rho_in = System.EquilibriumDensityMatrix(Hz_EnUnit,T,HT_approx)
rhoeq = rho_in.copy()
else:
    rho_in = np.sum(Sz,axis=0) # Initial Density Matrix
    rhoeq = np.sum(Sz,axis=0) # Equilibrium Density Matrix
    print("Trace of density metrix = ", np.trace(rho_in))


```

Trace of density metrix = 0j

In [74]: `Matrix(rho_in)`

Out[74]: $\begin{bmatrix} 0.5 & 0 \\ 0 & -0.5 \end{bmatrix}$

In [75]: `System.MatrixPlot(1,rho_in.real)`



```

/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source/PythonOnResonance.py:669: User
Warning: FixedFormatter should only be used together with FixedLocator
    ax.set_xticklabels([''] + labelx, fontsize=10)
/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source/PythonOnResonance.py:670: User
Warning: FixedFormatter should only be used together with FixedLocator
    ax.set_yticklabels([''] + labely, fontsize=10)

```

Zeeman Halitonian in Rotating Frame


In [76]: `OmegaRF = [-System.gammaH1*B0]`

```
Hzz = System.Zeeman_RotFrame(LarmorF, Sz, OmegaRF)
```

Pulse

```
In [77]: """  
Rotate the magnetization about Y-axis, by an angle theta.  
""";  
pulse_angle = 90.0  
rho = System.Rotate_H(rho_in,pulse_angle,np.sum(Sy,axis=0))
```

```
In [78]: System.MatrixPlot(2,rho.real)
```



Relaxation Constant

```
In [79]: '''  
Define longitudinal (R1) and transverse Relaxation (R2)  
R1 = [R1 of first spin, R1 of second spin,...]  
R2 = [R2 of first spin, R2 of second spin,...]  
''';  
  
R1 = np.asarray([0])  
R2 = np.asarray([0])  
System.Relaxation_Constants(R1,R2)  
  
'''
```

```
Options for "Rprocess": "No Relaxation" or "Phenomenological"
                        or "Random Field Fluxtuation" or "Dipolar"
'''
Rprocess = "No Relaxation"
```

Evolution of Density Matrix

```
In [80]: """
Sampling Rate, fs = n * Highest_Larmor_Frequency; minimum value of n = 2 (Nyquist-Shannon)
Dwell time, dt = 1/fs
Acquisition time, AQ is time for which we evolve the density matrix, in seconds.
Number of points in the simulation, Npoints
""";
Highest_Larmor_Frequency = 20.0
fs = 4 * Highest_Larmor_Frequency
dt = 1.0/fs
AQ = 5.0
Npoints = int(AQ/dt)
print("Number of points in the simulation", Npoints)

"""
option for solver, "method": "Unitary Propagator" or "ODE Solver"
"""
method = "Unitary Propagator"

start_time = time.time()
t, rho_t = System.Evolution_H(rhoeq, rho, Sx, Sy, Sz, Sp, Sm, Hzr, dt, Npoints, method, Rprocess)
end_time = time.time()
timetaken = end_time - start_time
print("Total time = %s seconds " % (timetaken))
```

Number of points in the simulation 400
Total time = 0.0039899349212646484 seconds

Expectation value

```
In [81]: """
Basis Operators in plus minus Z
""";
Basis = "PMZ spin half"
B_PMZ = System.SingleSpinOP(Sx, Sy, Sz, Sp, Sm, Basis)
```

Basis: $\frac{1}{\sqrt{2}}E, I_+, I_-, \sqrt{2}I_z$

```
In [82]: """
Lets see the expectation value of I+
""";

t, Mp = System.Expectation_H(rho_t, B_PMZ[1], dt, Npoints)
```

Windowing


```
In [83]: Mp = System.WindowFunction(t, Mp, 5.0)
```

Fourier Transform

```
In [84]: fs = 1.0/dt  
freq, spectrum = System.FourierTransform(Mp,fs,5)
```


Plotting

```
In [85]: System.Plotting(3,t,Mp,"time (s)","Mp","red")
```



```
/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part  
    return array(a, dtype, copy=False, order=order)  
No handles with labels found to put in legend.
```

```
In [86]: PH0 = 45.0  
spectrum_PH0 = System.PhaseAdjust_PH0(spectrum,PH0)  
System.Plotting(4,freq,spectrum_PH0,"Frequency (Hz)","Spectrum","blue")
```



```
/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part
  return array(a, dtype, copy=False, order=order)
No handles with labels found to put in legend.
```

Next tutorial: Evolution of Density Matrix in Hilbert Space Part 2

In this lecture you will see how to evolve the density matrix in time by solving Liouville-Von Neumann equation in Hilbert Space for two spin half system.

Any suggestion? write to me

If you see something is wrong please write to me, so that the PyOR can be error free.

vineethfrancis.physics@gmail.com

In []: