

Python On Resonance (PyOR)

Everybody can simulate NMR

Author: Vineeth Thalakkotloor

Email: vineethfrancis.physics@gmail.com

Tutorial 9: Incredible Natural Abundance Double Quantum Transfer Experiment (INADEQUATE)

INADEQUATE can be used to suppress solvent (H₂O) in H₁ spectroscopy. In this tutorial we consider three H₁ spins. Two are coupled and we will see how this pulse sequence will eliminate the third spin by phase cycling.

Some Basics:

Multiple Quantum Coherence → 2I_xS_x, 2I_yS_x, 2I_xS_y, 2I_yS_y

Double Quantum Coherence, DQ_x = 1/2 * (2I_xS_x - 2I_yS_y)

Double Quantum Coherence, DQ_y = 1/2 * (2I_yS_x + 2I_xS_y)

Zero Quantum Coherence, DQ_x = 1/2 * (2I_xS_x + 2I_yS_y)

Zero Quantum Coherence, DQ_y = 1/2 * (2I_yS_x - 2I_xS_y)

Reference book "NMR: The Toolkit, How Pulse Sequences Work" by P.J Hore, J.A. Jones and S. Wimperis

Load Python packages and define path to the source file "PythonOnResonance.py"

```
In [1]: pathSource = '/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source'
```

```
In [2]: from IPython.display import display, HTML
display(HTML("<style>.container { width:100% !important; }</style>"))
import sys
sys.path.append(pathSource)

import PythonOnResonance as PyOR

import time
import numpy as np
import matplotlib.pyplot as plt
from matplotlib import rc
%matplotlib notebook
import sympy as sp
from sympy import *
from IPython.display import display, Math, Latex
```

Generating Spin System

```
In [3]: """
Define Spin quantum numbers of your spins in "Slist1".
Slist1[0] is spin of first particle and Slist1[1] is spin of second particle.
""";

Slist1 = [1/2,1/2,1/2]
```

```
In [4]: """
Define Planck constant equals 1.
Because NMR spectroscopists are more interested to write Energy in frequency units.
if False then hbarEQ1 = hbar
""";

hbarEQ1 = True
```

```
In [5]: """
Generate Spin Operators
""";

System = PyOR.Numerical_MR(Slist1,hbarEQ1)

"""
Sx, Sy and Sz Operators
""";
Sx,Sy,Sz = System.SpinOperator()

"""
S+ and S- Operators
""";
Sp,Sm = System.PMoperators(Sx,Sy)
```

Zeeman Hamiltonian in Lab Frame

```
In [6]: """
Gyromagnetic Ratio
Gamma = [Gyromagnetic Ratio spin 1, Gyromagnetic Ratio spin 1, ...]
""";
Gamma = [System.gammaH1, System.gammaH1, System.gammaH1]

"""
Define the field of the spectrometer, B0 in Tesla.
""";
B0 = 9.4

"""
Define the chemical Shift of individual spins
Offset = [chemical Shift spin 1, chemical Shift spin 1, ..]
""";
Offset = [10,50,100] # Offset frequency in Hz

"""
Function "LarmorF" give the list Larmor frequencies of individual spins in lab frame
""";
LarmorF = System.LarmorFrequency(Gamma,B0,Offset)

Hz = System.Zeeman(LarmorF,Sz)
```

Larmor Frequency in MHz: [-400.22802765 -400.22806765 -400.22811765]

Initialize Density Matrix

In [7]:

```
"""
We will generate Initial Density Matrix in two ways:
First we will generate a density matrix as we prefer say, Sz.
Second we will create density matrix at thermal equilibrium

First Case
""";

Thermal_DensMatrix = False

if Thermal_DensMatrix:
    Hz_EnUnit = System.Convert_FreqUnitsTOEnergy(Hz)
    HT_approx = False # High Temperature Approximation is False
    T = 300 # Temperature in Kelvin
    rho_in = System.EquilibriumDensityMatrix(Hz_EnUnit,T,HT_approx)
    rhoeq = rho_in.copy()
else:
    rho_in = np.sum(Sz,axis=0) # Initial Density Matrix
    rhoeq = np.sum(Sz,axis=0) # Equilibrium Density Matrix
    print("Trace of density metrix = ", np.trace(rho_in))
```

Trace of density metrix = 0j

Zeeman Halitonian in Rotating Frame

In [8]:

```
OmegaRF = [-System.gammaH1*B0, -System.gammaH1*B0, -System.gammaH1*B0]
Hxr = System.Zeeman_RotFrame(LarmorF, Sz, OmegaRF)
```

J Coupling Hamiltonian

In [9]:

```
"""
Define J Coupling between each spins, Jlist[0][3] means J coupling between 1st spin and 4th spin
""";

Jlist = np.zeros((len(Slist1),len(Slist1)))
Jlist[0][1] = 10 # Spin 1 and 2 are coupled
Jlist[0][2] = 0 # spin 1 and 3 not coupled
Jlist[1][2] = 0 # spin 2 and 3 not coupled

Hj = System.Jcoupling_Weak(Jlist,Sz)
```

Relaxation Constant

In [10]:

```
'''
Define longitudinal (R1) and transverse Relaxation (R2)
''';

R1 = 1.0
R2 = 2.0
System.Relaxation_Constants(R1,R2)

'''
Options for "Rprocess": "No Relaxation" or "Phenomenological"
                        or "Random Field Fluxtuation" or "Dipolar"
''';
Rprocess = "No Relaxation"
```

Pulse 90x on all H1

```
In [11]: pulse_angle = 90.0
rho = System.Rotate_H(rho_in, pulse_angle, Sx[0]+Sx[1]+Sx[2])
```

Detect H1

```
In [12]: dt = 1.0e-4
fs = 1.0/dt
AQ = 5.0
Npoints = int(AQ/dt)
print("Number of points in the simulation", Npoints)

"""
option for solver, "method": "Unitary Propagator" or "ODE Solver"
"""
method = "Unitary Propagator"

start_time = time.time()
t, rho_t = System.Evolution_H(rhoeq, rho, Sx, Sy, Sz, Sp, Sm, Hzr + Hj, dt, Npoints, method, Rprocess)
end_time = time.time()
timetaken = end_time - start_time
print("Total time = %s seconds " % (timetaken))

Exp_H1 = (Sx[0]+Sx[1]+Sx[2]) + 1j * (Sy[0]+Sy[1]+Sy[2])
t, Mt_H1_normal = System.Expectation_H(rho_t, Exp_H1, dt, Npoints)
```

Number of points in the simulation 50000
Total time = 0.3216392993927002 seconds

INADEQUATE

```
In [13]: delta = 1/(4*Jlist[0][1])

# Phase of pulse for Phase Cycling
Ph1 = [0,90]
Ph2 = [90,180]
Ph3 = [0,90]
Sig_add = [1,-1]

def INADEQUATE(rho_in, rhoeq, delta, dt, AQ, Ph1, Ph2, Ph3, Sig_add):
    Npoints = int(AQ/dt)
    signal = np.zeros((Npoints), dtype=np.cdouble)
    method = "Unitary Propagator"
    l = 0
    for i,j,k in zip(Ph1, Ph2, Ph3):

        # 90 pulse on all spins
        rho = System.Rotate_H(rho_in, 90, System.Pulse_Phase(Sx, Sy, i))

        # First delta
        Dpoints = int(delta/dt)
        t, rho_t = System.Evolution_H(rhoeq, rho, Sx, Sy, Sz, Sp, Sm, Hzr + Hj, dt, Dpoints, method, Rprocess)

        # 180 pulse on all spins
        rho = System.Rotate_H(rho_t[-1], 180, System.Pulse_Phase(Sx, Sy, j))
```

```

# Second delta
t, rho_t = System.Evolution_H(rhoeq, rho, Sx, Sy, Sz, Sp, Sm, Hzr + Hj, dt, Dpoints, method,

# 90 pulse on all spins
rho = System.Rotate_H(rho_t[-1], 90, System.Pulse_Phase(Sx, Sy, k))

# 90 pulse on all spins
rho = System.Rotate_H(rho, 90, Sx[0]+Sx[1]+Sx[2])

# Detection
t, rho_t = System.Evolution_H(rhoeq, rho, Sx, Sy, Sz, Sp, Sm, Hzr + Hj, dt, Npoints, method,
Exp_H1 = (Sx[0]+Sx[1]+Sx[2]) + 1j * (Sy[0]+Sy[1]+Sy[2])
t, Mt_H1_inadequate = System.Expectation_H(rho_t, Exp_H1, dt, Npoints)

signal = signal + Sig_add[1] * Mt_H1_inadequate
l = l + 1

return t, signal

start_time = time.time()
t, Mt_H_inadequate = INADEQUATE(rho_in, rhoeq, delta, dt, AQ, Ph1, Ph2, Ph3, Sig_add)
end_time = time.time()
timetaken = end_time - start_time
print("Total time = %s seconds " % (timetaken))

```

Total time = 1.1896662712097168 seconds

Windowing

```

In [14]: Mt_H1_normal = System.WindowFunction(t, Mt_H1_normal, 5.0)
Mt_H_inadequate = System.WindowFunction(t, Mt_H_inadequate, 5.0)

```

Fourier Transform

```

In [15]: fs = 1.0/dt
freq, spectrum_normal = System.FourierTransform(Mt_H1_normal, fs, 5)
freq, spectrum_inadequate = System.FourierTransform(Mt_H_inadequate, fs, 5)


```

Plotting

```

In [16]: ...
Red: INADEQUATE
Blue: Normal (90deg pulse and detect H1)
...;
System.PlottingMulti(3, [t, t], [Mt_H_inadequate, Mt_H1_normal], "time (s)", R"M$_{T}$", ["red", '

```



```

/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part
  return array(a, dtype, copy=False, order=order)
/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part
  return array(a, dtype, copy=False, order=order)
No handles with labels found to put in legend.


```

In [17]:

```

...
Red: INADEQUATE
Blue: Normal (90deg pulse and detect H1)
''';
spectrum_inadequate_PH0 = System.PhaseAdjust_PH0(spectrum_inadequate,0.0)
Spectrum_normal_PH0 = System.PhaseAdjust_PH0(spectrum_normal,45.0)
System.PlottingMulti(4,[freq,freq],[spectrum_inadequate_PH0,Spectrum_normal_PH0],"Frequency

```



```

/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part
  return array(a, dtype, copy=False, order=order)

```

```
/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part
  return array(a, dtype, copy=False, order=order)
No handles with labels found to put in legend.
```

Any suggestion? write to me

If you see something is wrong please write to me, so that the PyOR can be error free.

vineethfrancis.physics@gmail.com

In []: