

# Python On Resonance (PyOR)

## Everybody can simulate NMR

Author: Vineeth Thalakkotloor

Email: vineethfrancis.physics@gmail.com

## Tutorial 5: Evolution of Density Matrix in Hilbert Space Part 2

In previous tutorial, we try to evolve the density matrix in time by solving Liouville-Von Neumann equation in Hilbert Space for a single spin half system. In this tutorial we will see two spin half system with out coupling.

### Load Python packages and define path to the source file "PythonOnResonance.py"

```
In [51]: pathSource = '/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source'
```

```
In [52]: from IPython.display import display, HTML
display(HTML("<style>.container { width:100% !important; }</style>"))
import sys
sys.path.append(pathSource)

import PythonOnResonance as PyOR

import time
import numpy as np
import matplotlib.pyplot as plt
from matplotlib import rc
%matplotlib notebook
import sympy as sp
from sympy import *
```

## Generating Spin System

```
In [53]: """
Define Spin quantum numbers of your spins in "Slist1".
Slist1[0] is spin of first particle and Slist1[1] is spin of second particle.
""";

Slist1 = [1/2,1/2]
```

```
In [54]: """
Define Planck constant equals 1.
Because NMR spectroscopists are more interested to write Energy in frequency units.
if False then hbarEQ1 = hbar
""";

hbarEQ1 = True
```

In [55]:

```
"""
Generate Spin Operators
""";

System = PyOR.Numerical_MR(Slist1,hbarEQ1)

"""
Sx, Sy and Sz Operators
""";
Sx,Sy,Sz = System.SpinOperator()

"""
S+ and S- Operators
""";
Sp,Sm = System.PMoperators(Sx,Sy)
```

## Zeeman Hamiltonian in Lab Frame

In [56]:

```
"""
Gyromagnetic Ratio
Gamma = [Gyromagnetic Ratio spin 1, Gyromagnetic Ratio spin 1, ...]
""";
Gamma = [System.gammaH1, System.gammaH1]

"""
Define the field of the spectrometer, B0 in Tesla.
""";
B0 = 9.4

"""
Define the chemical Shift of individual spins
Offset = [chemical Shift spin 1, chemical Shift spin 1, ..]
""";
Offset = [20,40] # Offset frequency in Hz

"""
Function "LarmorF" give the list Larmor frequencies of individual spins in lab frame
""";
LarmorF = System.LarmorFrequency(Gamma,B0,Offset)

Hz = System.Zeeman(LarmorF,Sz)
```

Larmor Frequency in MHz: [-400.22803765 -400.22805765]

## Initialize Density Matrix

In [57]:

```
"""
We will generate Initial Density Matrix in two ways:
First we will generate a density matrix as we prefer say, Sz.
Second we will create density matrix at thermal equilibrium

First Case
""";

Thermal_DensMatrix = False

if Thermal_DensMatrix:
    Hz_EnUnit = System.Convert_FreqUnitsTOEnergy(Hz)
    HT_approx = False # High Temperature Approximation is False
    T = 300 # Temperature in Kelvin
    rho_in = System.EquilibriumDensityMatrix(Hz_EnUnit,T,HT_approx)
```

```

rhoeq = rho_in.copy()
else:
    rho_in = np.sum(Sz,axis=0) # Initial Density Matrix
    rhoeq = np.sum(Sz,axis=0) # Equilibrium Density Matrix
    print("Trace of density metrix = ", np.trace(rho_in))

```


Trace of density metrix = 0j

In [58]: Matrix(rho\_in)

Out[58]:

$$\begin{bmatrix} 1.0 & 0 & 0 & 0 \\ 0 & 0 & 0 & 0 \\ 0 & 0 & 0 & 0 \\ 0 & 0 & 0 & -1.0 \end{bmatrix}$$

In [59]: System.MatrixPlot(1,rho\_in.real)



```

/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source/PythonOnResonance.py:669: User
Warning: FixedFormatter should only be used together with FixedLocator
    ax.set_xticklabels([''] + labelx,fontsize=10)
/media/HD2/Vineeth/PostDoc_Simulations/Github/PyOR_G/Source/PythonOnResonance.py:670: User
Warning: FixedFormatter should only be used together with FixedLocator
    ax.set_yticklabels([''] + labely,fontsize=10)

```


## Zeeman Halitonian in Rotating Frame

```
In [60]: OmegaRF = [-System.gammaH1*B0, -System.gammaH1*B0]
Hrz = System.Zeeman_RotFrame(LarmorF, Sz, OmegaRF)
```

## Pulse

```
In [61]: """
Rotate the magnetization about Y-axis, by an angle theta.
""";
pulse_angle = 90.0
rho = System.Rotate_H(rho_in, pulse_angle, np.sum(Sy, axis=0))
```

```
In [62]: System.MatrixPlot(2, rho.real)
```



## Relaxation Constant

```
In [63]: ...
Define longitudinal (R1) and transverse Relaxation (R2)
R1 = [R1 of first spin, R1 of second spin,...]
R2 = [R2 of first spin, R2 of second spin,...]
''';

R1 = np.asarray([0,0])
R2 = np.asarray([0,0])
System.Relaxation_Constants(R1,R2)
```

```

'''
Options for "Rprocess": "No Relaxation" or "Phenomenological"
                        or "Random Field Fluxtuation" or "Dipolar"
'''
Rprocess = "No Relaxation"

```

## Evolution of Density Matrix

In [64]:

```

"""
Sampling Rate, fs = n * Highest_Larmor_Frequency; minimum value of n = 2 (Nyquist-Shannon)
Dwell time, dt = 1/fs
Acquisition time, AQ is time for which we evolve the density matrix, in seconds.
Number of points in the simulation, Npoints
"""
Highest_Larmor_Frequency = 40.0
fs = 4 * Highest_Larmor_Frequency
dt = 1.0/fs
AQ = 5.0
Npoints = int(AQ/dt)
print("Number of points in the simulation", Npoints)

"""
option for solver, "method": "Unitary Propagator" or "ODE Solver"
"""
method = "Unitary Propagator"

start_time = time.time()
t, rho_t = System.Evolution_H(rhoeq, rho, Sx, Sy, Sz, Sp, Sm, Hzr, dt, Npoints, method, Rprocess)
end_time = time.time()
timetaken = end_time - start_time
print("Total time = %s seconds " % (timetaken))

```

Number of points in the simulation 800  
Total time = 0.010149955749511719 seconds

## Expectation value

In [65]:

```

"""
Basis Operators in plus minus Z
"""
Basis = "PMZ spin half"
B_PMZ = System.TwoSpinOP(Sx, Sy, Sz, Sp, Sm, Basis)

```

Basis:  $\frac{1}{2}E, \frac{1}{\sqrt{2}}I_+, \frac{1}{\sqrt{2}}I_-, \frac{1}{\sqrt{2}}I_z, \frac{1}{\sqrt{2}}S_+, \frac{1}{\sqrt{2}}S_-, S_z, \sqrt{2}I_xS_z, \sqrt{2}I_yS_z, \sqrt{2}I_zS_x, \sqrt{2}I_zS_y, 2I_zS_z, I_+S_+, I_+S_-, I_-S_+, I_-S_-$

In [66]:

```

"""
Lets see the expectation value of I+ and S+
"""

t, Mp = System.Expectation_H(rho_t, B_PMZ[1] + B_PMZ[4], dt, Npoints)

```

## Windowing

In [67]:

```

Mp = System.WindowFunction(t, Mp, 5.0)


```

# Fourier Transform

```
In [68]: fs = 1.0/dt  
         freq, spectrum = System.FourierTransform(Mp,fs,5)
```


## Plotting

```
In [69]: System.Plotting(3,t,Mp,"time (s)","Mp","red")
```



```
/opt/anaconda3/lib/python3.9/site-packages/numpy/core/_asarray.py:102: ComplexWarning: Casting complex values to real discards the imaginary part  
    return array(a, dtype, copy=False, order=order)  
No handles with labels found to put in legend.
```

```
In [71]: System.Plotting(4,freq,np.absolute(spectrum),"Frequency (Hz)","Spectrum","blue")
```



No handles with labels found to put in legend.

## Next tutorial: Evolution of Density Matrix in Hilbert Space Part 3

In this lecture you will see how to evolve the density matrix in time by solving Liouville-Von Neumann equation in Hilbert Space of 6 spin half system (ethanol).

### Any suggestion? write to me

If you see something is wrong please write to me, so that the PyOR can be error free.

[vineethfrancis.physics@gmail.com](mailto:vineethfrancis.physics@gmail.com)

In [ ]: